

## Product datasheet for **TS417480P5**

### Laminin beta 2 (LAMB2) CytoSection

#### Product data:

Product Type:	CytoSections
Description:	Transient overexpression of LAMB2 in HEK293T cells, FFPE control for IHC, ICC and ISH staining, 25 slides per pack
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC217480
Tag:	C-MYC/DDK
Detection Antibodies:	Clone OTI4C5, Anti-DDK (FLAG) monoclonal antibody (TA50011-100)
Target Detection Antibodies:	Laminin beta 2 (LAMB2) Mouse Monoclonal Antibody [Clone ID: OTI3B4] (TA808897)
ACCN:	<a href="#">NM_002292</a> , <a href="#">NP_002283</a>
Synonyms:	LAMS; NPHS5
Storage:	Room Temperature
Stability:	Slides are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin. 5 µm sections of the FFPE cell pellet blocks are cut and mounted on positively charged SuperFrost slides.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	<a href="#">NP_002283</a>
Locus ID:	3913
Cytogenetics:	3p21.31
Protein Families:	Druggable Genome, Secreted Protein
Protein Pathways:	ECM-receptor interaction, Focal adhesion, Pathways in cancer, Small cell lung cancer



[View online »](#)