

Product datasheet for **TS417522P5**

TRIF (TICAM1) CytoSection

Product data:

Product Type:	CytoSections
Description:	Transient overexpression of TICAM1 in HEK293T cells, FFPE control for IHC, ICC and ISH staining, 25 slides per pack
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC217522
Tag:	C-MYC/DDK
Detection Antibodies:	DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)
Target Detection Antibodies:	TRIF (TICAM1) Mouse Monoclonal Antibody [Clone ID: OT11G7] (TA800679)
ACCN:	NM_182919 , NP_891549
Synonyms:	IIAE6; MyD88-3; PRVTIRB; TICAM-1; TRIF
Storage:	Room Temperature
Stability:	Slides are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin. 5 µm sections of the FFPE cell pellet blocks are cut and mounted on positively charged SuperFrost slides.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	NP_891549
Locus ID:	148022
Cytogenetics:	19p13.3
Protein Families:	Druggable Genome
Protein Pathways:	Toll-like receptor signaling pathway



[View online »](#)