

## Product datasheet for **TS417559P5**

### Glutamyl Prolyl tRNA synthetase (EPRS) CytoSection

#### Product data:

Product Type:	CytoSections
Description:	Transient overexpression of EPRS1 in HEK293T cells, FFPE control for IHC, ICC and ISH staining, 25 slides per pack
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC217559
Tag:	C-MYC/DDK
Detection Antibodies:	Clone OTI4C5, Anti-DDK (FLAG) monoclonal antibody (TA50011-100)
ACCN:	<a href="#">NM_004446</a> , <a href="#">NP_004437</a>
Synonyms:	EARS; EPRS; GLUPRORS; HLD15; PARS; PIG32; QARS; QPRS
Storage:	Room Temperature
Stability:	Slides are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin. 5 µm sections of the FFPE cell pellet blocks are cut and mounted on positively charged SuperFrost slides.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	<a href="#">NP_004437</a>
Locus ID:	2058
Cytogenetics:	1q41
Protein Families:	Druggable Genome
Protein Pathways:	Aminoacyl-tRNA biosynthesis, Metabolic pathways, Porphyrin and chlorophyll metabolism



[View online »](#)