

## Product datasheet for **TS417752P5**

### DC2L1 (DYNC2LI1) CytoSection

#### Product data:

Product Type:	CytoSections
Description:	Transient overexpression of DYNC2LI1, transcript variant 3, in HEK293T cells, FFPE control for IHC, ICC and ISH staining, 25 slides per pack
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC217752
Tag:	C-MYC/DDK
Detection Antibodies:	DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)
ACCN:	<u>NM_001012665</u> , <u>NP_001012683</u>
Synonyms:	CGI-60; D2LIC; DKFZp564A033; dynein, cytoplasmic 2, light intermediate chain 1; dynein 2 light intermediate chain; LIC3
Storage:	Room Temperature
Stability:	Slides are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin. 5 µm sections of the FFPE cell pellet blocks are cut and mounted on positively charged SuperFrost slides.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	<u>NP_001012683</u>
Locus ID:	51626
Cytogenetics:	2p21



[View online »](#)