

Product datasheet for **TS417827P5**

ADAM28 CytoSection

Product data:

Product Type:	CytoSections
Description:	Transient overexpression of ADAM28, transcript variant 3, in HEK293T cells, FFPE control for IHC, ICC and ISH staining, 25 slides per pack
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC217827
Tag:	C-MYC/DDK
Detection Antibodies:	DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)
Target Detection Antibodies:	ADAM28 Mouse Monoclonal Antibody [Clone ID: OTI3F4] (TA812739)
ACCN:	<u>NM_021777</u> , <u>NP_068547</u>
Synonyms:	ADAM 28; eMDC II; eMDCII; MDC-L; MDCL
Storage:	Room Temperature
Stability:	Slides are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin. 5 µm sections of the FFPE cell pellet blocks are cut and mounted on positively charged SuperFrost slides.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	<u>NP_068547</u>
Locus ID:	10863
Cytogenetics:	8p21.2
Protein Families:	Druggable Genome, Protease, Secreted Protein, Transmembrane



[View online »](#)