

Product datasheet for **TS417979P5**

ASAP1 CytoSection

Product data:

Product Type:	CytoSections
Description:	Transient overexpression of ASAP1 in HEK293T cells, FFPE control for IHC, ICC and ISH staining, 25 slides per pack
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC217979
Tag:	C-MYC/DDK
Detection Antibodies:	Clone OTI4C5, Anti-DDK (FLAG) monoclonal antibody (TA50011-100)
ACCN:	NM_018482 , NP_060952
Synonyms:	AMAP1; CENTB4; DDEF1; PAG2; PAP; ZG14P
Storage:	Room Temperature
Stability:	Slides are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin. 5 µm sections of the FFPE cell pellet blocks are cut and mounted on positively charged SuperFrost slides.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	NP_060952
Locus ID:	50807
Cytogenetics:	8q24.21-q24.22
Protein Pathways:	Endocytosis, Fc gamma R-mediated phagocytosis



[View online »](#)