

Product datasheet for **TS418082**

GH2 CytoSection

Product data:

| | |
|--|---|
| Product Type: | CytoSections |
| Description: | Transient overexpression of GH2, transcript variant 4, in HEK293T cells, FFPE control for IHC, ICC and ISH staining, 5 slides per pack |
| Species: | Human |
| Expression Host: | HEK293T |
| Expression cDNA Clone or AA Sequence: | TrueORF Clone RC218082 |
| Tag: | C-MYC/DDK |
| Detection Antibodies: | DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569) |
| Target Detection Antibodies: | GH2 Mouse Monoclonal Antibody [Clone ID: OT15B11] (TA811356) |
| ACCN: | <u>NM_022556</u> , <u>NP_072050</u> |
| Synonyms: | GH-V; GHB2; GHL; GHV; hGH-V |
| Storage: | Room Temperature |
| Stability: | Slides are guaranteed for a year from the date of receipt if proper storage instructions were followed. |
| Preparation: | HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin. 5 µm sections of the FFPE cell pellet blocks are cut and mounted on positively charged SuperFrost slides. |
| Note: | This product is for research use only and is not approved for use in humans or in clinical diagnosis. |
| RefSeq: | <u>NP_072050</u> |
| Locus ID: | 2689 |
| Cytogenetics: | 17q23.3 |
| Protein Families: | Druggable Genome, Secreted Protein |
| Protein Pathways: | Cytokine-cytokine receptor interaction, Jak-STAT signaling pathway, Neuroactive ligand-receptor interaction |



[View online »](#)