

Product datasheet for **TS418101P5**

MEK3 (MAP2K3) CytoSection

Product data:

| | |
|---------------------------------------|---|
| Product Type: | CytoSections |
| Description: | Transient overexpression of MAP2K3, transcript variant A, in HEK293T cells, FFPE control for IHC, ICC and ISH staining, 25 slides per pack |
| Species: | Human |
| Expression Host: | HEK293T |
| Expression cDNA Clone or AA Sequence: | TrueORF Clone RC218101 |
| Tag: | C-MYC/DDK |
| Detection Antibodies: | DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569) |
| Target Detection Antibodies: | MEK3 (MAP2K3) Mouse Monoclonal Antibody [Clone ID: OTI4B3] (TA505701) |
| ACCN: | NM_002756 , NP_002747 |
| Synonyms: | MAPKK3; MEK3; MKK3; PRKMK3; SAPKK-2; SAPKK2 |
| Storage: | Room Temperature |
| Stability: | Slides are guaranteed for a year from the date of receipt if proper storage instructions were followed. |
| Preparation: | HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin. 5 µm sections of the FFPE cell pellet blocks are cut and mounted on positively charged SuperFrost slides. |
| Note: | This product is for research use only and is not approved for use in humans or in clinical diagnosis. |
| RefSeq: | NP_002747 |
| Locus ID: | 5606 |
| Cytogenetics: | 17p11.2 |
| Protein Families: | Druggable Genome, Protein Kinase, Transcription Factors |
| Protein Pathways: | Amyotrophic lateral sclerosis (ALS), Fc epsilon RI signaling pathway, GnRH signaling pathway, MAPK signaling pathway, Toll-like receptor signaling pathway |



[View online »](#)