

## Product datasheet for **TS418141**

### TGF alpha (TGFA) CytoSection

#### Product data:

Product Type:	CytoSections
Description:	Transient overexpression of TGFA, transcript variant 2, in HEK293T cells, FFPE control for IHC, ICC and ISH staining, 5 slides per pack
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC218141
Tag:	C-MYC/DDK
Detection Antibodies:	DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)
ACCN:	<u><a href="#">NM_001099691</a></u> , <u><a href="#">NP_001093161</a></u>
Synonyms:	TFGA
Storage:	Room Temperature
Stability:	Slides are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin. 5 µm sections of the FFPE cell pellet blocks are cut and mounted on positively charged SuperFrost slides.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	<u><a href="#">NP_001093161</a></u>
Locus ID:	7039
Cytogenetics:	2p13.3
Protein Families:	Druggable Genome, Secreted Protein, Transmembrane
Protein Pathways:	ErbB signaling pathway, Glioma, Non-small cell lung cancer, Pancreatic cancer, Pathways in cancer, Prostate cancer, Renal cell carcinoma



[View online »](#)