

## **Product datasheet for TS418162P5**

## OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

## **CSPS (SULT1A3) CytoSection**

**Product data:** 

**Product Type:** CytoSections

**Description:** Transient overexpression of SULT1A3 in HEK293T cells, FFPE control for IHC, ICC and ISH

staining, 25 slides per pack

Species: Human
Expression Host: HEK293T

Expression cDNA Clone

or AA Sequence:

TrueORF Clone RC218162

Tag: C-MYC/DDK

**Detection Antibodies:** Clone OTI4C5, Anti-DDK (FLAG) monoclonal antibody (TA50011-100)

**Target Detection** 

**Antibodies:** 

CSPS (SULT1A3) Mouse Monoclonal Antibody [Clone ID: OTI3G5] (TA811313)

ACCN: <u>NM 177552</u>, <u>NP 808220</u>

Synonyms: HAST; HAST3; M-PST; ST1A3; ST1A3/ST1A4; ST1A4; ST1A5; STM; TL-PST

**Storage:** Room Temperature

**Stability:** Slides are guaranteed for a year from the date of receipt if proper storage instructions were

followed.

**Preparation:** HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells

were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin &

dehydrated before embedding in paraffin. 5 µm sections of the FFPE cell pellet blocks are cut

and mounted on positively charged SuperFrost slides.

**Note:** This product is for research use only and is not approved for use in humans or in clinical

diagnosis.

RefSeq: NP 808220

**Locus ID:** 6818

Cytogenetics: 16p11.2

**Protein Pathways:** Sulfur metabolism

