

## Product datasheet for **TS418310P5**


### PSMA (FOLH1) CytoSection

#### Product data:

Product Type:	CytoSections
Description:	Transient overexpression of FOLH1, transcript variant 1, in HEK293T cells, FFPE control for IHC, ICC and ISH staining, 25 slides per pack
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC218310
Tag:	C-MYC/DDK
Detection Antibodies:	DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)
Target Detection Antibodies:	PSMA (FOLH1) Mouse Monoclonal Antibody [Clone ID: OTI3H5] (TA504752)
ACCN:	<a href="#">NM_004476</a> , <a href="#">NP_004467</a>
Synonyms:	FGCP; FOLH; GCP2; GCPII; mGCP; NAALAD1; NAALAdase; PSM; PSMA
Storage:	Room Temperature
Stability:	Slides are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin. 5 µm sections of the FFPE cell pellet blocks are cut and mounted on positively charged SuperFrost slides.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	<a href="#">NP_004467</a>
Locus ID:	2346
Cytogenetics:	11p11.12
Protein Families:	Druggable Genome, Protease, Transmembrane



[View online »](#)

**Product images:**

Immunohistochemistry staining on FOLH1 overexpressed cytosection TS418310P5 with rabbit anti DDK clone OTIR5G2 C/N [TA592569] at 1:2000 dilution 20m RT. Antibody staining was achieved with HIER Citrate pH6 , Polink1 with DAB chromogen detection kit ([D11-18]). Positive stain shown with the brown chromogen present. HEK293T cells were transfected with cDNA clone [RC218310], 5 micron sections, 40x magnification