

## Product datasheet for TS418323P5

## OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

## **ALAS2 CytoSection**

**Product data:** 

**Product Type:** CytoSections

Description: Transient overexpression of ALAS2, transcript variant 3, in HEK293T cells, FFPE control for

IHC, ICC and ISH staining, 25 slides per pack

Species: Human **Expression Host:** HEK293T

**Expression cDNA Clone** 

or AA Sequence:

TrueORF Clone RC218323

C-MYC/DDK Tag:

**Detection Antibodies:** DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)

ALAS2 Mouse Monoclonal Antibody [Clone ID: OTI2D3] (TA809959)

**Target Detection** 

Antibodies:

ACCN:

NM 001037968, NP 001033057

ALAS-E; ALASE; ANH1; ASB; SIDBA1; XLDPP; XLEPP; XLSA Synonyms:

Storage: Room Temperature

Stability: Slides are guaranteed for a year from the date of receipt if proper storage instructions were

followed.

Preparation: HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells

were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin &

dehydrated before embedding in paraffin. 5 µm sections of the FFPE cell pellet blocks are cut

and mounted on positively charged SuperFrost slides.

This product is for research use only and is not approved for use in humans or in clinical Note:

diagnosis.

RefSeq: NP 001033057

Locus ID: 212

Cytogenetics: Xp11.21

**Protein Families:** Druggable Genome

**Protein Pathways:** Glycine, serine and threonine metabolism, Metabolic pathways, Porphyrin and chlorophyll

metabolism

