

Product datasheet for TS418490P5

CD1C CytoSection

Product data:

Product Type: CytoSections **Description:** Transient overexpression of CD1C in HEK293T cells, FFPE control for IHC, ICC and ISH staining, 25 slides per pack Species: Human **Expression Host:** HEK293T **Expression cDNA Clone** TrueORF Clone RC218490 or AA Sequence: C-MYC/DDK Tag: **Detection Antibodies:** DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569) **Target Detection** CD1C Mouse Monoclonal Antibody [Clone ID: OTI2F4] (TA505411) Antibodies: ACCN: NM 001765, NP 001756 BDCA1; CD1; CD1A; R7 Synonyms: Storage: **Room Temperature** Stability: Slides are guaranteed for a year from the date of receipt if proper storage instructions were followed. **Preparation:** HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin. 5 µm sections of the FFPE cell pellet blocks are cut and mounted on positively charged SuperFrost slides. This product is for research use only and is not approved for use in humans or in clinical Note: diagnosis. **RefSeq:** NP 001756 Locus ID: 911 **Cytogenetics:** 1q23.1 **Protein Families:** Druggable Genome, Transmembrane **Protein Pathways:** Hematopoietic cell lineage




This product is to be used for laboratory only. Not for diagnostic or therapeutic use. ©2023 OriGene Technologies, Inc., 9620 Medical Center Drive, Ste 200, Rockville, MD 20850, US

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

Product images:



Immunohistology staining on CD1C overexpressed cytosection TS418490P5 with mouse anti DDK clone OTI11C3 C/N [TA180144] at 1:400 dilution 4C O/N. Antibody staining was achieved with HIER Citrate pH6 , Polink1 with DAB chromogen detection kit ([D11-18]). Positive stain shown with the brown chromogen present. HEK293T cells were transfected with cDNA clone [RC218490], 5 micron sections, 40x magnification

This product is to be used for laboratory only. Not for diagnostic or therapeutic use. ©2023 OriGene Technologies, Inc., 9620 Medical Center Drive, Ste 200, Rockville, MD 20850, US