

Product datasheet for **TS418544**

Myosin Heavy chain 1 (MYH1) CytoSection

Product data:

Product Type:	CytoSections
Description:	Transient overexpression of MYH1 in HEK293T cells, FFPE control for IHC, ICC and ISH staining, 25 slides per pack
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC218544
Tag:	C-MYC/DDK
Detection Antibodies:	DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)
ACCN:	NM_005963 , NP_005954
Synonyms:	HEL71; MYHa; MyHC-2x; MyHC-2X/D; MYHSA1
Storage:	Room Temperature
Stability:	Slides are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin. 5 µm sections of the FFPE cell pellet blocks are cut and mounted on positively charged SuperFrost slides.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	NP_005954
Locus ID:	4619
Cytogenetics:	17p13.1
Protein Pathways:	Tight junction, Viral myocarditis



[View online »](#)