

## Product datasheet for **TS419188**

### VPS13B CytoSection

#### Product data:

Product Type:	CytoSections
Description:	Transient overexpression of VPS13B, transcript variant 4, in HEK293T cells, FFPE control for IHC, ICC and ISH staining, 5 slides per pack
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC219188
Tag:	C-MYC/DDK
Detection Antibodies:	DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)
Target Detection Antibodies:	VPS13B Mouse Monoclonal Antibody [Clone ID: OTI6F6] (TA807801)
ACCN:	<a href="#">NM_181661</a> , <a href="#">NP_858047</a>
Synonyms:	CHS1; COH1
Storage:	Room Temperature
Stability:	Slides are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin. 5 µm sections of the FFPE cell pellet blocks are cut and mounted on positively charged SuperFrost slides.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	<a href="#">NP_858047</a>
Locus ID:	157680
Cytogenetics:	8q22.2



[View online »](#)