

## Product datasheet for **TS419398P5**

### Angiopoietin 1 (ANGPT1) CytoSection

#### Product data:

Product Type:	CytoSections
Description:	Transient overexpression of ANGPT1, transcript variant 1, in HEK293T cells, FFPE control for IHC, ICC and ISH staining, 25 slides per pack
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC219398
Tag:	C-MYC/DDK
Detection Antibodies:	DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)
Target Detection Antibodies:	Angiopoietin 1 (ANGPT1) Mouse Monoclonal Antibody [Clone ID: OT17F9] (TA804638)
ACCN:	<u><a href="#">NM_001146</a></u> , <u><a href="#">NP_001137</a></u>
Synonyms:	AGP1; AGPT; ANG1; HAE5
Storage:	Room Temperature
Stability:	Slides are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin. 5 µm sections of the FFPE cell pellet blocks are cut and mounted on positively charged SuperFrost slides.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	<u><a href="#">NP_001137</a></u>
Locus ID:	284
Cytogenetics:	8q23.1
Protein Families:	Druggable Genome, ES Cell Differentiation/IPS, Secreted Protein



[View online »](#)