

## Product datasheet for **TS419524**

### MS4A7 CytoSection

#### Product data:

Product Type:	CytoSections
Description:	Transient overexpression of MS4A7, transcript variant 4, in HEK293T cells, FFPE control for IHC, ICC and ISH staining, 5 slides per pack
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC219524
Tag:	C-MYC/DDK
Detection Antibodies:	DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)
ACCN:	<a href="#">NM_206940</a> , <a href="#">NP_996823</a>
Synonyms:	4SPAN2; CD20L4; CFFM4; MS4A8
Storage:	Room Temperature
Stability:	Slides are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin. 5 $\mu$ m sections of the FFPE cell pellet blocks are cut and mounted on positively charged SuperFrost slides.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	<a href="#">NP_996823</a>
Locus ID:	58475
Cytogenetics:	11q12.2
Protein Families:	Druggable Genome, Transmembrane



[View online »](#)