

Product datasheet for **TS419579P5**

NIRF (UHRF2) CytoSection

Product data:

Product Type:	CytoSections
Description:	Transient overexpression of UHRF2 in HEK293T cells, FFPE control for IHC, ICC and ISH staining, 25 slides per pack
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC219579
Tag:	C-MYC/DDK
Detection Antibodies:	DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)
Target Detection Antibodies:	NIRF (UHRF2) Mouse Monoclonal Antibody [Clone ID: OTI4E10] (TA810202)
ACCN:	NM_152896 , NP_690856
Synonyms:	NIRF; RNF107; TDRD23; URF2
Storage:	Room Temperature
Stability:	Slides are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin. 5 µm sections of the FFPE cell pellet blocks are cut and mounted on positively charged SuperFrost slides.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	NP_690856
Locus ID:	115426
Cytogenetics:	9p24.1
Protein Families:	Druggable Genome, Transcription Factors



[View online »](#)