

Product datasheet for **TS419668P5**

Aquaporin 8 (AQP8) CytoSection

Product data:

Product Type:	CytoSections
Description:	Transient overexpression of AQP8 in HEK293T cells, FFPE control for IHC, ICC and ISH staining, 25 slides per pack
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC219668
Tag:	C-MYC/DDK
Detection Antibodies:	Clone OTI4C5, Anti-DDK (FLAG) monoclonal antibody (TA50011-100)
Target Detection Antibodies:	Aquaporin 8 (AQP8) Mouse Monoclonal Antibody [Clone ID: OTI6C9] (TA810776)
ACCN:	<u>NM_001169</u> , <u>NP_001160</u>
Synonyms:	AQP-8
Storage:	Room Temperature
Stability:	Slides are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin. 5 µm sections of the FFPE cell pellet blocks are cut and mounted on positively charged SuperFrost slides.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	<u>NP_001160</u>
Locus ID:	343
Cytogenetics:	16p12.1
Protein Families:	Druggable Genome, Transmembrane



[View online »](#)