

## Product datasheet for **TS419723**

### CAMK2B CytoSection

#### Product data:

Product Type:	CytoSections
Description:	Transient overexpression of CAMK2B, transcript variant 4, in HEK293T cells, FFPE control for IHC, ICC and ISH staining, 5 slides per pack
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC219723
Tag:	C-MYC/DDK
Detection Antibodies:	DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)
Target Detection Antibodies:	CAMK2B Mouse Monoclonal Antibody [Clone ID: OT18D9] (TA808345)
ACCN:	<u><a href="#">NM_172080</a></u> , <u><a href="#">NP_742077</a></u>
Synonyms:	CAM2; CAMK2; CAMKB; CaMKIIbeta; MRD54
Storage:	Room Temperature
Stability:	Slides are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin. 5 µm sections of the FFPE cell pellet blocks are cut and mounted on positively charged SuperFrost slides.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	<u><a href="#">NP_742077</a></u>
Locus ID:	816
Cytogenetics:	7p13
Protein Families:	Druggable Genome, Protein Kinase



[View online »](#)

**Protein Pathways:**

Calcium signaling pathway, ErbB signaling pathway, Glioma, GnRH signaling pathway, Long-term potentiation, Melanogenesis, Neurotrophin signaling pathway, Olfactory transduction, Oocyte meiosis, Wnt signaling pathway