

## Product datasheet for **TS419808**

### DR5 (TNFRSF10B) CytoSection

#### Product data:

Product Type:	CytoSections
Description:	Transient overexpression of TNFRSF10B, transcript variant 2, in HEK293T cells, FFPE control for IHC, ICC and ISH staining, 5 slides per pack
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC219808
Tag:	C-MYC/DDK
Detection Antibodies:	DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)
Target Detection Antibodies:	DR5 (TNFRSF10B) Mouse Monoclonal Antibody [Clone ID: OT11B9] (TA807279)
ACCN:	<u><a href="#">NM_147187</a></u> , <u><a href="#">NP_671716</a></u>
Synonyms:	CD262; DR5; KILLER; KILLER/DR5; TRAIL-R2; TRAILR2; TRICK2; TRICK2A; TRICK2B; TRICKB; ZTNFR9
Storage:	Room Temperature
Stability:	Slides are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin. 5 µm sections of the FFPE cell pellet blocks are cut and mounted on positively charged SuperFrost slides.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	<u><a href="#">NP_671716</a></u>
Locus ID:	8795
Cytogenetics:	8p21.3
Protein Families:	Druggable Genome, Transmembrane



[View online »](#)

**Protein Pathways:** Apoptosis, Cytokine-cytokine receptor interaction, Natural killer cell mediated cytotoxicity, p53 signaling pathway