

Product datasheet for **TS419842P5**

SCNN1D CytoSection

Product data:

Product Type:	CytoSections
Description:	Transient overexpression of SCNN1D, transcript variant 2, in HEK293T cells, FFPE control for IHC, ICC and ISH staining, 25 slides per pack
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC219842
Tag:	C-MYC/DDK
Detection Antibodies:	DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)
ACCN:	NM_002978 , NP_002969
Synonyms:	dNaCh; ENaCd; ENaCdelta; MGC149710; MGC149711; OTTHUMP00000002084; SCNED; sodium channel, nonvoltage-gated 1, delta; sodium channel, voltage-gated, type I, delta polypeptide
Storage:	Room Temperature
Stability:	Slides are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin. 5 µm sections of the FFPE cell pellet blocks are cut and mounted on positively charged SuperFrost slides.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	NP_002969
Locus ID:	6339
Cytogenetics:	1p36.33
Protein Families:	Druggable Genome, Ion Channels: Other, Transmembrane



[View online »](#)