

## Product datasheet for **TS420235**

### Glucocorticoid Receptor (NR3C1) CytoSection

#### Product data:

Product Type:	CytoSections
Description:	Transient overexpression of NR3C1, transcript variant 4, in HEK293T cells, FFPE control for IHC, ICC and ISH staining, 5 slides per pack
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC220235
Tag:	C-MYC/DDK
Detection Antibodies:	DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)
Target Detection Antibodies:	Glucocorticoid Receptor (NR3C1) Mouse Monoclonal Antibody [Clone ID: OTI6B5] (TA806259)
ACCN:	<a href="#">NM_001018076</a> , <a href="#">NP_001018086</a>
Synonyms:	GCCR; GCR; GCRST; GR; GRL
Storage:	Room Temperature
Stability:	Slides are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin. 5 µm sections of the FFPE cell pellet blocks are cut and mounted on positively charged SuperFrost slides.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	<a href="#">NP_001018086</a>
Locus ID:	2908
Cytogenetics:	5q31.3
Protein Families:	Druggable Genome, Nuclear Hormone Receptor, Transcription Factors
Protein Pathways:	Neuroactive ligand-receptor interaction



[View online »](#)