

Product datasheet for **TS420387P5**

eRF1 (ETF1) CytoSection

Product data:

Product Type:	CytoSections
Description:	Transient overexpression of ETF1 in HEK293T cells, FFPE control for IHC, ICC and ISH staining, 25 slides per pack
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC220387
Tag:	C-MYC/DDK
Detection Antibodies:	DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)
Target Detection Antibodies:	eRF1 (ETF1) Mouse Monoclonal Antibody [Clone ID: OTI2E3] (TA502891)
ACCN:	<u>NM_004730</u> , <u>NP_004721</u>
Synonyms:	D5S1995; ERF; ERF1; RF1; SUP45L1; TB3-1
Storage:	Room Temperature
Stability:	Slides are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin. 5 µm sections of the FFPE cell pellet blocks are cut and mounted on positively charged SuperFrost slides.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	<u>NP_004721</u>
Locus ID:	2107
Cytogenetics:	5q31.2


[View online »](#)