

Product datasheet for **TS420513P5**

HNRPM (HNRNPM) CytoSection

Product data:

Product Type:	CytoSections
Description:	Transient overexpression of HNRNPM, transcript variant 2, in HEK293T cells, FFPE control for IHC, ICC and ISH staining, 25 slides per pack
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC220513
Tag:	C-MYC/DDK
Detection Antibodies:	DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)
Target Detection Antibodies:	HNRPM (HNRNPM) Mouse Monoclonal Antibody [Clone ID: OT13F3] (TA803107)
ACCN:	<u>NM_031203</u> , <u>NP_112480</u>
Synonyms:	CEAR; hnRNP M; HNRNPM4; HNRPM; HNRPM4; HTGR1; NAGR1
Storage:	Room Temperature
Stability:	Slides are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin. 5 µm sections of the FFPE cell pellet blocks are cut and mounted on positively charged SuperFrost slides.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	<u>NP_112480</u>
Locus ID:	4670
Cytogenetics:	19p13.2
Protein Families:	Druggable Genome
Protein Pathways:	Spliceosome



[View online »](#)