

Product datasheet for **TS420572P5**

AK3L1 (AK4) CytoSection

Product data:

Product Type:	CytoSections
Description:	Transient overexpression of AK4, transcript variant 6, in HEK293T cells, FFPE control for IHC, ICC and ISH staining, 25 slides per pack
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC220572
Tag:	C-MYC/DDK
Detection Antibodies:	DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)
Target Detection Antibodies:	AK3L1 (AK4) Mouse Monoclonal Antibody [Clone ID: OTI3B1] (TA503011)
ACCN:	NM_013410 , NP_037542
Synonyms:	AK3; AK3L1; AK3L2; AK 4
Storage:	Room Temperature
Stability:	Slides are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin. 5 µm sections of the FFPE cell pellet blocks are cut and mounted on positively charged SuperFrost slides.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	NP_037542
Locus ID:	205
Cytogenetics:	1p31.3
Protein Families:	Druggable Genome
Protein Pathways:	Metabolic pathways, Purine metabolism



[View online »](#)