

## Product datasheet for **TS420582P5**

### ASIC3 CytoSection

#### Product data:

Product Type:	CytoSections
Description:	Transient overexpression of ASIC3, transcript variant 1, in HEK293T cells, FFPE control for IHC, ICC and ISH staining, 25 slides per pack
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC220582
Tag:	C-MYC/DDK
Detection Antibodies:	DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)
ACCN:	<u>NM_004769</u> , <u>NP_004760</u>
Synonyms:	ACCN3; DRASIC; SLNAC1; TNaC1
Storage:	Room Temperature
Stability:	Slides are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin. 5 µm sections of the FFPE cell pellet blocks are cut and mounted on positively charged SuperFrost slides.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	<u>NP_004760</u>
Locus ID:	9311
Cytogenetics:	7q36.1
Protein Families:	Druggable Genome, Ion Channels: Other



[View online »](#)