

## Product datasheet for **TS420827P5**

### ABCG8 CytoSection

#### Product data:

|                                       |                                                                                                                                                                                                                                                                                                                                        |
|---------------------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Type:                         | CytoSections                                                                                                                                                                                                                                                                                                                           |
| Description:                          | Transient overexpression of ABCG8 in HEK293T cells, FFPE control for IHC, ICC and ISH staining, 25 slides per pack                                                                                                                                                                                                                     |
| Species:                              | Human                                                                                                                                                                                                                                                                                                                                  |
| Expression Host:                      | HEK293T                                                                                                                                                                                                                                                                                                                                |
| Expression cDNA Clone or AA Sequence: | TrueORF Clone RC220827                                                                                                                                                                                                                                                                                                                 |
| Tag:                                  | C-MYC/DDK                                                                                                                                                                                                                                                                                                                              |
| Detection Antibodies:                 | Clone OTI4C5, Anti-DDK (FLAG) monoclonal antibody (TA50011-100)                                                                                                                                                                                                                                                                        |
| ACCN:                                 | <u><a href="#">NM_022437</a></u> , <u><a href="#">NP_071882</a></u>                                                                                                                                                                                                                                                                    |
| Synonyms:                             | GBD4; STSL; STSL1                                                                                                                                                                                                                                                                                                                      |
| Storage:                              | Room Temperature                                                                                                                                                                                                                                                                                                                       |
| Stability:                            | Slides are guaranteed for a year from the date of receipt if proper storage instructions were followed.                                                                                                                                                                                                                                |
| Preparation:                          | HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin. 5 $\mu$ m sections of the FFPE cell pellet blocks are cut and mounted on positively charged SuperFrost slides. |
| Note:                                 | This product is for research use only and is not approved for use in humans or in clinical diagnosis.                                                                                                                                                                                                                                  |
| RefSeq:                               | <u><a href="#">NP_071882</a></u>                                                                                                                                                                                                                                                                                                       |
| Locus ID:                             | 64241                                                                                                                                                                                                                                                                                                                                  |
| Cytogenetics:                         | 2p21                                                                                                                                                                                                                                                                                                                                   |
| Protein Families:                     | Druggable Genome, Transmembrane                                                                                                                                                                                                                                                                                                        |
| Protein Pathways:                     | ABC transporters                                                                                                                                                                                                                                                                                                                       |



[View online »](#)