

## Product datasheet for TS420914P5

## OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

## **MBD1 CytoSection**

**Product data:** 

**Product Type:** CytoSections

Description: Transient overexpression of MBD1, transcript variant 1, in HEK293T cells, FFPE control for

IHC, ICC and ISH staining, 25 slides per pack

Species: Human **Expression Host:** HEK293T

**Expression cDNA Clone** 

or AA Sequence:

TrueORF Clone RC220914

C-MYC/DDK Tag:

**Detection Antibodies:** DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)

**Target Detection** 

Antibodies:

MBD1 Mouse Monoclonal Antibody [Clone ID: OTI2H5] (TA811767)

ACCN: NM 015846, NP 056671

CXXC3; PCM1; RFT Synonyms: Storage: Room Temperature

Stability: Slides are guaranteed for a year from the date of receipt if proper storage instructions were

followed.

Preparation: HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells

were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin &

dehydrated before embedding in paraffin. 5 µm sections of the FFPE cell pellet blocks are cut

and mounted on positively charged SuperFrost slides.

This product is for research use only and is not approved for use in humans or in clinical Note:

diagnosis.

RefSeq: NP 056671

Locus ID: 4152

Cytogenetics: 18q21.1

**Protein Families:** Druggable Genome, Transcription Factors

