

Product datasheet for **TS421010**

Kallikrein 8 (KLK8) CytoSection

Product data:

Product Type:	CytoSections
Description:	Transient overexpression of KLK8, transcript variant 2, in HEK293T cells, FFPE control for IHC, ICC and ISH staining, 5 slides per pack
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC221010
Tag:	C-MYC/DDK
Detection Antibodies:	DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)
Target Detection Antibodies:	Kallikrein 8 (KLK8) Mouse Monoclonal Antibody [Clone ID: OT11A9] (TA506137)
ACCN:	NM_144505 , NP_653088
Synonyms:	HNP; NP; NRPN; PRSS19; TADG14
Storage:	Room Temperature
Stability:	Slides are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin. 5 µm sections of the FFPE cell pellet blocks are cut and mounted on positively charged SuperFrost slides.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	NP_653088
Locus ID:	11202
Cytogenetics:	19q13.41
Protein Families:	Druggable Genome, Secreted Protein, Transmembrane



[View online »](#)