

Product datasheet for **TS421018P5**

NUMB CytoSection

Product data:

Product Type:	CytoSections
Description:	Transient overexpression of NUMB, transcript variant 2, in HEK293T cells, FFPE control for IHC, ICC and ISH staining, 25 slides per pack
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC221018
Tag:	C-MYC/DDK
Detection Antibodies:	DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)
Target Detection Antibodies:	NUMB Mouse Monoclonal Antibody [Clone ID: OT11A7] (TA501614)
ACCN:	NM_001005744 , NP_001005744
Synonyms:	C14orf41; c14_5527; S171
Storage:	Room Temperature
Stability:	Slides are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin. 5 µm sections of the FFPE cell pellet blocks are cut and mounted on positively charged SuperFrost slides.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	NP_001005744
Locus ID:	8650
Cytogenetics:	14q24.2-q24.3
Protein Pathways:	Notch signaling pathway



[View online »](#)