

Product datasheet for **TS421140P5**

Caveolin 3 (CAV3) CytoSection

Product data:

Product Type:	CytoSections
Description:	Transient overexpression of CAV3, transcript variant 1, in HEK293T cells, FFPE control for IHC, ICC and ISH staining, 25 slides per pack
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC221140
Tag:	C-MYC/DDK
Detection Antibodies:	DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)
ACCN:	NM_033337 , NP_203123
Synonyms:	LGMD1C; LQT9; MPDT; RMD2; VIP-21; VIP21
Storage:	Room Temperature
Stability:	Slides are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin. 5 µm sections of the FFPE cell pellet blocks are cut and mounted on positively charged SuperFrost slides.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	NP_203123
Locus ID:	859
Cytogenetics:	3p25.3
Protein Families:	Druggable Genome, Transmembrane
Protein Pathways:	Focal adhesion



[View online »](#)