

## OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

## Product datasheet for TS421303P5

## **Progesterone Receptor (PGR) CytoSection**

## **Product data:**

| Product Type:                            | CytoSections  |
|--|---|
| Description:                             | Transient overexpression of PGR, transcript variant 2, in HEK293T cells, FFPE control for IHC,<br>ICC and ISH staining, 25 slides per pack  |
| Species:                                 | Human   |
| Expression Host:                         | HEK293T   |
| Expression cDNA Clone<br>or AA Sequence: | TrueORF Clone RC221303  |
| Tag:                                     | C-MYC/DDK   |
| <b>Detection Antibodies:</b>             | DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)  |
| Target Detection<br>Antibodies:          | Progesterone Receptor (PGR) Mouse Monoclonal Antibody [Clone ID: OTI11E8] (TA802704)  |
| ACCN:                                    | <u>NM 000926, NP 000917</u>   |
| Synonyms:                                | NR3C3; PR   |
| Storage:                                 | Room Temperature  |
| Stability:                               | Slides are guaranteed for a year from the date of receipt if proper storage instructions were followed.   |
| Preparation:                             | HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin. 5 μm sections of the FFPE cell pellet blocks are cut and mounted on positively charged SuperFrost slides. |
| Note:                                    | This product is for research use only and is not approved for use in humans or in clinical diagnosis.   |
| RefSeq:                                  | <u>NP 000917</u>  |
| Locus ID:                                | 5241  |
| Cytogenetics:                            | 11q22.1   |
| Protein Families:                        | Druggable Genome, Nuclear Hormone Receptor, Transcription Factors   |
| Protein Pathways:                        | Oocyte meiosis, Progesterone-mediated oocyte maturation   |



This product is to be used for laboratory only. Not for diagnostic or therapeutic use. ©2023 OriGene Technologies, Inc., 9620 Medical Center Drive, Ste 200, Rockville, MD 20850, US