

Product datasheet for **TS421319P5**

C3orf15 (MAATS1) CytoSection

Product data:

| | |
|--|---|
| Product Type: | CytoSections |
| Description: | Transient overexpression of CFAP91 in HEK293T cells, FFPE control for IHC, ICC and ISH staining, 25 slides per pack |
| Species: | Human |
| Expression Host: | HEK293T |
| Expression cDNA Clone or AA Sequence: | TrueORF Clone RC221319 |
| Tag: | C-MYC/DDK |
| Detection Antibodies: | DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569) |
| Target Detection Antibodies: | C3orf15 (MAATS1) Mouse Monoclonal Antibody [Clone ID: OTI6G1] (TA802361) |
| ACCN: | <u>NM_033364</u> , <u>NP_203528</u> |
| Synonyms: | AAT1; AAT1alpha; C3orf15; CaM-IP2; MAATS1; SPATA26 |
| Storage: | Room Temperature |
| Stability: | Slides are guaranteed for a year from the date of receipt if proper storage instructions were followed. |
| Preparation: | HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin. 5 µm sections of the FFPE cell pellet blocks are cut and mounted on positively charged SuperFrost slides. |
| Note: | This product is for research use only and is not approved for use in humans or in clinical diagnosis. |
| RefSeq: | <u>NP_203528</u> |
| Locus ID: | 89876 |
| Cytogenetics: | 3q13.33 |



[View online »](#)