

## Product datasheet for **TS421341**

### AGAP1 CytoSection

#### Product data:

Product Type:	CytoSections
Description:	Transient overexpression of AGAP1, transcript variant 1, in HEK293T cells, FFPE control for IHC, ICC and ISH staining, 5 slides per pack
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC221341
Tag:	C-MYC/DDK
Detection Antibodies:	DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)
Target Detection Antibodies:	AGAP1 Mouse Monoclonal Antibody [Clone ID: OT14E8] (TA504204)
ACCN:	<a href="#">NM_001037131</a> , <a href="#">NP_001032208</a>
Synonyms:	AGAP-1; CENTG2; cnt-g2; GGAP1
Storage:	Room Temperature
Stability:	Slides are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin. 5 µm sections of the FFPE cell pellet blocks are cut and mounted on positively charged SuperFrost slides.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	<a href="#">NP_001032208</a>
Locus ID:	116987
Cytogenetics:	2q37.2
Protein Pathways:	Endocytosis



[View online »](#)