

Product datasheet for **TS421478P5**

Junctional Adhesion Molecule 1 (F11R) CytoSection

Product data:

Product Type:	CytoSections
Description:	Transient overexpression of F11R in HEK293T cells, FFPE control for IHC, ICC and ISH staining, 25 slides per pack
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC221478
Tag:	C-MYC/DDK
Detection Antibodies:	Clone OTI4C5, Anti-DDK (FLAG) monoclonal antibody (TA50011-100)
Target Detection Antibodies:	Junctional Adhesion Molecule 1 (F11R) Mouse Monoclonal Antibody [Clone ID: OTI3H3] (TA506017)
ACCN:	NM_016946 , NP_058642
Synonyms:	CD321; JAM; JAM1; JAMA; JCAM; KAT; PAM-1
Storage:	Room Temperature
Stability:	Slides are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin. 5 µm sections of the FFPE cell pellet blocks are cut and mounted on positively charged SuperFrost slides.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	NP_058642
Locus ID:	50848
Cytogenetics:	1q23.3
Protein Families:	Druggable Genome, Transmembrane
Protein Pathways:	Cell adhesion molecules (CAMs), Epithelial cell signaling in Helicobacter pylori infection, Leukocyte transendothelial migration, Tight junction



[View online »](#)