

Product datasheet for **TS421499P5**

SOAT 2 (SOAT2) CytoSection

Product data:

Product Type:	CytoSections
Description:	Transient overexpression of SOAT2 in HEK293T cells, FFPE control for IHC, ICC and ISH staining, 25 slides per pack
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC221499
Tag:	C-MYC/DDK
Detection Antibodies:	Clone OTI4C5, Anti-DDK (FLAG) monoclonal antibody (TA50011-100)
Target Detection Antibodies:	SOAT 2 (SOAT2) Mouse Monoclonal Antibody [Clone ID: OTI2D3] (TA800775)
ACCN:	<u>NM_003578</u> , <u>NP_003569</u>
Synonyms:	ACACT2; ACAT2; ARGP2
Storage:	Room Temperature
Stability:	Slides are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin. 5 µm sections of the FFPE cell pellet blocks are cut and mounted on positively charged SuperFrost slides.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	<u>NP_003569</u>
Locus ID:	8435
Cytogenetics:	12q13.13
Protein Families:	Transmembrane
Protein Pathways:	Steroid biosynthesis



[View online »](#)