

Product datasheet for **TS421637P5**

MAP2 CytoSection

Product data:

Product Type:	CytoSections
Description:	Transient overexpression of MAP2, transcript variant 2, in HEK293T cells, FFPE control for IHC, ICC and ISH staining, 25 slides per pack
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC221637
Tag:	C-MYC/DDK
Detection Antibodies:	DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)
Target Detection Antibodies:	MAP2 Mouse Monoclonal Antibody [Clone ID: OT15D1] (TA813550)
ACCN:	NM_031845 , NP_114033
Synonyms:	MAP-2; MAP2A; MAP2B; MAP2C
Storage:	Room Temperature
Stability:	Slides are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin. 5 µm sections of the FFPE cell pellet blocks are cut and mounted on positively charged SuperFrost slides.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	NP_114033
Locus ID:	4133
Cytogenetics:	2q34
Protein Families:	Adult stem cells, Druggable Genome, ES Cell Differentiation/IPS



[View online »](#)