

## Product datasheet for **TS421707P5**


### PCDHA11 CytoSection

#### Product data:

Product Type:	CytoSections
Description:	Transient overexpression of PCDHA11, transcript variant 1, in HEK293T cells, FFPE control for IHC, ICC and ISH staining, 25 slides per pack
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC221707
Tag:	C-MYC/DDK
Detection Antibodies:	DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)
ACCN:	<u>NM_018902</u> , <u>NP_061725</u>
Synonyms:	CNR7; CNRN7; CNRS7; CRNR7; PCDH-ALPHA11
Storage:	Room Temperature
Stability:	Slides are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin. 5 µm sections of the FFPE cell pellet blocks are cut and mounted on positively charged SuperFrost slides.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	<u>NP_061725</u>
Locus ID:	56138
Cytogenetics:	5q31.3
Protein Families:	Transmembrane



[View online »](#)

**Product images:**

Immunohistology staining on PCDHA11 overexpressed cytosection TS421707P5 with rabbit anti DDK clone OTIR5G2 C/N [TA592569] at 1:2000 dilution 20m RT. Antibody staining was achieved with HIER Citrate pH6 , Polink1 with DAB chromogen detection kit ([D11-18]). Positive stain shown with the brown chromogen present. HEK293T cells were transfected with cDNA clone [RC221707], 5 micron sections, 40x magnification