

Product datasheet for **TS421838P5**

TrkB (NTRK2) CytoSection

Product data:

Product Type:	CytoSections
Description:	Transient overexpression of NTRK2, transcript variant c, in HEK293T cells, FFPE control for IHC, ICC and ISH staining, 25 slides per pack
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC221838
Tag:	C-MYC/DDK
Detection Antibodies:	DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)
Target Detection Antibodies:	TrkB (NTRK2) Mouse Monoclonal Antibody [Clone ID: OTI2E1] (TA500386)
ACCN:	NM_001018064 , NP_001018074
Synonyms:	DEE58; EIEE58; GP145-TrkB; OBHD; trk-B; TRKB
Storage:	Room Temperature
Stability:	Slides are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin. 5 µm sections of the FFPE cell pellet blocks are cut and mounted on positively charged SuperFrost slides.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	NP_001018074
Locus ID:	4915
Cytogenetics:	9q21.33
Protein Families:	Druggable Genome, Protein Kinase, Transmembrane
Protein Pathways:	MAPK signaling pathway, Neurotrophin signaling pathway



[View online »](#)