

Product datasheet for TS421860P5

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

IL1 Receptor II (IL1R2) CytoSection

Product data:

Product Type: CytoSections

Description: Transient overexpression of IL1R2, transcript variant 1, in HEK293T cells, FFPE control for IHC,

ICC and ISH staining, 25 slides per pack

Species: Human
Expression Host: HEK293T

Expression cDNA Clone

or AA Sequence:

TrueORF Clone RC221860

Tag: C-MYC/DDK

Detection Antibodies: DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)

Target Detection

Antibodies:

IL1 Receptor II (IL1R2) Mouse Monoclonal Antibody [Clone ID: OTI3C6] (TA506802)

ACCN: NM 004633, NP 004624

Synonyms: CD121b; CDw121b; IL-1R-2; IL-1RT-2; IL-1RT2; IL1R2c; IL1RB

Storage: Room Temperature

Stability: Slides are guaranteed for a year from the date of receipt if proper storage instructions were

followed.

Preparation: HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells

were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin &

dehydrated before embedding in paraffin. 5 µm sections of the FFPE cell pellet blocks are cut

and mounted on positively charged SuperFrost slides.

Note: This product is for research use only and is not approved for use in humans or in clinical

diagnosis.

RefSeq: <u>NP 004624</u>

Locus ID: 7850 **Cytogenetics:** 2q11.2

Protein Families: Druggable Genome, Transmembrane

Protein Pathways: Cytokine-cytokine receptor interaction, Hematopoietic cell lineage, MAPK signaling pathway

