

Product datasheet for **TS421884**

CK1 epsilon (CSNK1E) CytoSection

Product data:

Product Type:	CytoSections
Description:	Transient overexpression of CSNK1E, transcript variant 2, in HEK293T cells, FFPE control for IHC, ICC and ISH staining, 5 slides per pack
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC221884
Tag:	C-MYC/DDK
Detection Antibodies:	DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)
Target Detection Antibodies:	CK1 epsilon (CSNK1E) Mouse Monoclonal Antibody [Clone ID: OTI5D4] (TA803799)
ACCN:	NM_001894 , NP_001885
Synonyms:	CK1e; CK1epsilon; HCK1E
Storage:	Room Temperature
Stability:	Slides are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin. 5 µm sections of the FFPE cell pellet blocks are cut and mounted on positively charged SuperFrost slides.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	NP_001885
Locus ID:	1454
Cytogenetics:	22q13.1
Protein Families:	Druggable Genome, Protein Kinase
Protein Pathways:	Circadian rhythm - mammal, Hedgehog signaling pathway, Wnt signaling pathway



[View online »](#)