

Product datasheet for **TS422464P5**

NEUROD4 CytoSection

Product data:

Product Type:	CytoSections
Description:	Transient overexpression of NEUROD4 in HEK293T cells, FFPE control for IHC, ICC and ISH staining, 25 slides per pack
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC222464
Tag:	C-MYC/DDK
Detection Antibodies:	Clone OTI4C5, Anti-DDK (FLAG) monoclonal antibody (TA50011-100)
Target Detection Antibodies:	NEUROD4 Mouse Monoclonal Antibody [Clone ID: OTI1B7] (TA810658)
ACCN:	NM_021191 , NP_067014
Synonyms:	ATH-3; ATH3; Atoh3; bHLHa4; MATH-3; MATH3
Storage:	Room Temperature
Stability:	Slides are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin. 5 µm sections of the FFPE cell pellet blocks are cut and mounted on positively charged SuperFrost slides.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	NP_067014
Locus ID:	58158
Cytogenetics:	12q13.2
Protein Families:	Druggable Genome, Transcription Factors



[View online »](#)