

Product datasheet for **TS422669P5**

SLC19A1 CytoSection

Product data:

Product Type:	CytoSections
Description:	Transient overexpression of SLC19A1, transcript variant 1, in HEK293T cells, FFPE control for IHC, ICC and ISH staining, 25 slides per pack
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC222669
Tag:	C-MYC/DDK
Detection Antibodies:	DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)
ACCN:	<u>NM_194255</u> , <u>NP_919231</u>
Synonyms:	CHMD; FOLT; hRFC; hSLC19A1; IFC-1; IFC1; MEGAF; REFC; RFC; RFC1; RFT-1
Storage:	Room Temperature
Stability:	Slides are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin. 5 µm sections of the FFPE cell pellet blocks are cut and mounted on positively charged SuperFrost slides.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	<u>NP_919231</u>
Locus ID:	6573
Cytogenetics:	21q22.3
Protein Families:	Stem cell - Pluripotency, Transmembrane



[View online »](#)