

Product datasheet for **TS422716**

MTGR1 (CBFA2T2) CytoSection

Product data:

Product Type:	CytoSections
Description:	Transient overexpression of CBFA2T2, transcript variant 2, in HEK293T cells, FFPE control for IHC, ICC and ISH staining, 5 slides per pack
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC222716
Tag:	C-MYC/DDK
Detection Antibodies:	DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)
ACCN:	NM_005093 , NP_005084
Synonyms:	EHT; MTGR1; p85; ZMYND3
Storage:	Room Temperature
Stability:	Slides are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin. 5 μ m sections of the FFPE cell pellet blocks are cut and mounted on positively charged SuperFrost slides.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	NP_005084
Locus ID:	9139
Cytogenetics:	20q11.21-q11.22
Protein Families:	Transcription Factors



[View online »](#)