

Product datasheet for **TS422840P5**

CYP27B1 CytoSection

Product data:

Product Type:	CytoSections
Description:	Transient overexpression of CYP27B1 in HEK293T cells, FFPE control for IHC, ICC and ISH staining, 25 slides per pack
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC222840
Tag:	C-MYC/DDK
Detection Antibodies:	Clone OTI4C5, Anti-DDK (FLAG) monoclonal antibody (TA50011-100)
ACCN:	NM_000785 , NP_000776
Synonyms:	CP2B; CYP1; CYP1alpha; CYP27B; P450c1; PDDR; VDD1; VDDR; VDDRI; VDR
Storage:	Room Temperature
Stability:	Slides are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin. 5 µm sections of the FFPE cell pellet blocks are cut and mounted on positively charged SuperFrost slides.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	NP_000776
Locus ID:	1594
Cytogenetics:	12q14.1
Protein Families:	Druggable Genome, P450
Protein Pathways:	Metabolic pathways, Steroid biosynthesis



[View online »](#)