

Product datasheet for **TS422951**

Aprataxin (APTX) CytoSection

Product data:

Product Type:	CytoSections
Description:	Transient overexpression of APTX, transcript variant 4, in HEK293T cells, FFPE control for IHC, ICC and ISH staining, 5 slides per pack
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC222951
Tag:	C-MYC/DDK
Detection Antibodies:	DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)
ACCN:	NM_017692 , NP_060162
Synonyms:	AOA; AOA1; aprataxin; ataxia 1, early onset with hypoalbuminemia; AXA1; EAOH; EOAHA; FHA-HIT; FLJ20157; MGC1072; OTTHUMP00000021188
Storage:	Room Temperature
Stability:	Slides are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin. 5 µm sections of the FFPE cell pellet blocks are cut and mounted on positively charged SuperFrost slides.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	NP_060162
Locus ID:	54840
Cytogenetics:	9p21.1
Protein Families:	Druggable Genome



[View online »](#)