

## Product datasheet for **TS422951P5**

### Aprataxin (APTX) CytoSection

#### Product data:

|                                       |                                                                                                                                                                                                                                                                                                                                   |
|---------------------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Type:                         | CytoSections                                                                                                                                                                                                                                                                                                                      |
| Description:                          | Transient overexpression of APTX, transcript variant 4, in HEK293T cells, FFPE control for IHC, ICC and ISH staining, 25 slides per pack                                                                                                                                                                                          |
| Species:                              | Human                                                                                                                                                                                                                                                                                                                             |
| Expression Host:                      | HEK293T                                                                                                                                                                                                                                                                                                                           |
| Expression cDNA Clone or AA Sequence: | TrueORF Clone RC222951                                                                                                                                                                                                                                                                                                            |
| Tag:                                  | C-MYC/DDK                                                                                                                                                                                                                                                                                                                         |
| Detection Antibodies:                 | DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)                                                                                                                                                                                                                                                |
| ACCN:                                 | <u><a href="#">NM_017692</a></u> , <u><a href="#">NP_060162</a></u>                                                                                                                                                                                                                                                               |
| Synonyms:                             | AOA; AOA1; aprataxin; ataxia 1, early onset with hypoalbuminemia; AXA1; EAOH; EOAHA; FHA-HIT; FLJ20157; MGC1072; OTTHUMP00000021188                                                                                                                                                                                               |
| Storage:                              | Room Temperature                                                                                                                                                                                                                                                                                                                  |
| Stability:                            | Slides are guaranteed for a year from the date of receipt if proper storage instructions were followed.                                                                                                                                                                                                                           |
| Preparation:                          | HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin. 5 µm sections of the FFPE cell pellet blocks are cut and mounted on positively charged SuperFrost slides. |
| Note:                                 | This product is for research use only and is not approved for use in humans or in clinical diagnosis.                                                                                                                                                                                                                             |
| RefSeq:                               | <u><a href="#">NP_060162</a></u>                                                                                                                                                                                                                                                                                                  |
| Locus ID:                             | 54840                                                                                                                                                                                                                                                                                                                             |
| Cytogenetics:                         | 9p21.1                                                                                                                                                                                                                                                                                                                            |
| Protein Families:                     | Druggable Genome                                                                                                                                                                                                                                                                                                                  |



[View online »](#)