

Product datasheet for **TS423133P5**

UBE2V1 (TMEM189-UBE2V1) CytoSection

Product data:

Product Type:	CytoSections
Description:	Transient overexpression of PEDS1-UBE2V1 in HEK293T cells, FFPE control for IHC, ICC and ISH staining, 25 slides per pack
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC223133
Tag:	C-MYC/DDK
Detection Antibodies:	DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)
ACCN:	NM_199203 , NP_954673
Synonyms:	CROC-1B; CROC1B; KUA-UEV; TMEM189-UBE2V1
Storage:	Room Temperature
Stability:	Slides are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin. 5 μ m sections of the FFPE cell pellet blocks are cut and mounted on positively charged SuperFrost slides.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	NP_954673
Locus ID:	387522
Cytogenetics:	20q13.13
Protein Families:	Transmembrane



[View online »](#)