

## Product datasheet for **TS423136P5**

### NMDAR2A (GRIN2A) CytoSection

#### Product data:

Product Type:	CytoSections
Description:	Transient overexpression of GRIN2A, transcript variant 2, in HEK293T cells, FFPE control for IHC, ICC and ISH staining, 25 slides per pack
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC223136
Tag:	C-MYC/DDK
Detection Antibodies:	DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)
ACCN:	<a href="#">NM_000833</a> , <a href="#">NP_000824</a>
Synonyms:	EPND; FESD; GluN2A; LKS; NMDAR2A; NR2A
Storage:	Room Temperature
Stability:	Slides are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin. 5 µm sections of the FFPE cell pellet blocks are cut and mounted on positively charged SuperFrost slides.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	<a href="#">NP_000824</a>
Locus ID:	2903
Cytogenetics:	16p13.2
Protein Families:	Druggable Genome, Ion Channels: Glutamate Receptors, Ion Channels: Sodium, Transmembrane
Protein Pathways:	Alzheimer's disease, Amyotrophic lateral sclerosis (ALS), Calcium signaling pathway, Long-term potentiation, Neuroactive ligand-receptor interaction, Systemic lupus erythematosus



[View online »](#)